LEVAGE

Reshaping Laser System

New Generation Laser Endolifting



FacialContouring™

It is an outpatient procedure for facial and body remodeling, with tightening and volume reconfiguration effects, treating signs of aging without pain and on an outpatient basis.

Why Levage?

- Minimally invasive or without anesthesia
- Instant results.
- A single treatment with long-term effects.
- Non-surgical
- Combinable with aesthetic procedures.





√ REMODELS VOLUME

ANTES / DESPUÉS















Unique Advantages



High ResistanceHighly durable aluminum.



WIFI UpdatesFree update to the latest



5 Years Warranty
On the diode and 3 years the unit.



Instructional VideosMemory for video files.



Why 1470nm?

The 1470nm wavelength offers a greater affinity for fat and water than other lasers, therefore emulsifying fat faster and with less energy. This results in minimal recovery time, painlessness and lower risk of damage than other laser wavelengths.

Accessories







Tips



Lipo Kit



Fiber Aspiration



Kit 200,400,600











Especificaciones

Wavelength
Maximum power

Pulse width

Screen/Software

Fiber Connection
Operation mode

Weight

Dimensiones

1470 nm

12W

10-25000 ms/1-25 s

Pantalla táctil vívida de 7 HD/sistema

operativo inteligente

200 con conector estándar SMA905

Modo continuo/Modo de pulso

2KG NW.

218 nm x 47 nm x 155 nm